



Recognizing, Evaluating, and Minimizing Artifacts

Robertson Chapter 10

Introduction

- Webster's Medical definition = ECG and EEG wave that arises from sources other than heart or brain
 - Too basic of a definition for PSG
- Early recognition and evaluation is essential
 - Must be corrected in timely and accurate fashion
 - Also must determine when artifacts can't be eliminated without jeopardizing data
- Artifact-free recording starts with proper electrode prep and placement



Introduction

- Ensure accurate signal processing through using a quality system and optimal amplifier settings
- Vigilant watch of the recording integrity is important



Electrode and Sensor Application



- Typically, the tech applies:
 - EEG, EOG, ECG and EMG electrodes
 - RIP belts
 - Airflow pressure sensor
 - Thermal airflow sensor
 - Snore sensor
 - Oximeter probe

Electrode Application Methods



- Collodion
 - Pros:
 - Very stable application
 - Fewer artifacts
 - Cons:
 - Highly flammable
 - Irritating to upper respiratory tract and mucous membranes
 - Should be applied in well ventilated area
 - Must be stored in cool, dry area
 - Often in a vented fireproof container

Electrode Application Methods

- Conductive paste
 - Pros:
 - Nontoxic
 - Requires no special storage or air-handling procedures
 - Cons:
 - Not as stable as collodion
 - Increased likelihood for artifacts
 - Use gauze to hold head electrodes in place
 - Use hypoallergenic tape for face and body



Electrode Removal

- Acetone is used to remove collodion
 - Must be stored in cool, dry location
 - May need to be kept in vented, fireproof cabinet
 - Be careful around patient's eyes, nose, and mouth
 - Be careful if patient has respiratory disease because of the fumes
- Paste removal requires no special chemicals
 - Best for patient to take a hot shower after study, let hot water run over them for a few minutes before applying shampoo

Gloving

- Gloves should always be used during direct patient contact
 - Use during electrode site prep, electrode application, and electrode removal
- If latex gloves are used, check with patient regarding allergies
 - Should have non-latex alternatives
- May have to wear multiple layers of gloves when using acetone and collodion as they can break down gloves
- Use new gloves for each patient

Electrode Site Preparation

- Assess patient's skin
 - Is the patient sensitive to adhesives?
 - Does the patient have a skin condition?
- Prepare skin using skin prep lotion with abrasive material
 - Helps remove oils, skin care products, and dead skin cells
 - Only prep contact area for electrode
 - Prepping larger area can lead to signal contamination and salt bridge artifact
 - Salt bridging = Excessive conductive material and electrodes in close proximity communicating
 - Be careful not to be too excessive with abrasion

Electrode Composition

- Electrodes should be composed of gold, silver, or silver-silver chloride
- Gold cup electrodes are gold-plated silver
 - Be careful not to scratch or crack as this causes artifact
- Avoid mixing electrode materials within a derivation
 - Will lead to distortion of recorded data
 - Leads to electrical potential differences

Electrode Impedance

- Check prior to and during acquisition
- Should be under 5,000 Ohms
- If impedances are imbalanced, CMR becomes ineffective
- Low impedance value
 - Very low = May be contact with another electrode due to excess gel, paste, or skin prep, sweat, or another electrolyte liquid
 - Check if impedance less than 1,000 Ohms

Electrode Impedance

- High impedance value
 - Need more conductive material between electrode and skin
 - Electrode needs cleaned
 - Bad electrode wire
 - Electrode popping
 - Makes recording more prone to 60 Hz interference
- Make sure connections are secure in head box

Differential Amplifier

- Amplifies difference between 2 input signals
 - Electrical difference between its 2 inputs
- Typically housed in head box
- CMRR = Ability of an amplifier to reject in-phase potentials and amplify out-of-phase potentials
 - Defining feature of differential amplifier
 - Input 1 = Exploring electrode
 - Input 2 = Reference electrode

Filters

- Allow manipulation of signals to attenuate artifacts
 - Should only be used as a final option as they change appearance of data
- LFF
 - Attenuates frequencies approaching, at, and below cut-off frequency
 - If at 1 Hz, a 1 Hz input will show 80% of original amplitude
 - Faster signals will show > 80% but < 100% of original amplitude
 - Less than 1 Hz will show < 80% of original amplitude
 - Changes time for amplitude to decay back to baseline


Filters

- HFF
 - Attenuates frequencies approaching, at, or above cut-off frequency
 - Effect mirrors LFF
- 50/60 Hz Filter
 - Designed to minimize data associated with electrical current interference
 - Also known as “line filter” or “notch filter” as it is designed to “notch out” frequencies surrounding line noise

Artifact Isolation

- Determine if issue is in more than 1 channel
- Is the artifact mechanical or physiological?
- Is the electrode loose?
- Is the lead wire intact?
- Is there enough conductive material?
- Is there a short circuit?

What is the Source of the Artifact?




Line Frequency Artifact

- Electrical current noise
- Common artifact
 - May be unavoidable if current is strong
 - Usually seen in bioelectric channels but sometimes seen in transduced channels
 - If seen in all bioelectric channels, typically is a bad ground electrode
- May appear in 1 channel or all channels
 - Uniform in appearance
- Activating 60 Hz filter will make this artifact disappear but should be last resort

Line Frequency Artifact


- If impedance good, amplifier may be picking up ambient leakage from nearby sources
 - Typically this affects all channels
- AASM recommends that line filter not be used on EMG
- Not all fast artifacts are 60-Hz line noise

Ambient Signals with Head Box Disconnected



Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

High Frequency Artifact on RAT




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Muscle Artifact


- Typically appears in EEG, EOG, and occasionally ECG
- Generated by localized muscle activity/ tension
 - Has frequency of 20-200 Hz
 - May look a lot like 60 Hz artifact
 - May disappear at sleep onset
- Abdominal effort
 - Reduce gain/sensitivity so signals are not overlapping

Respiratory Event with Muscle Artifact




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Leg Cramps



Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Leg Cramps




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Bruxism




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

EEG Arousal



Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Body Rocking




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

ECG Artifact

- Common artifact
- Typically seen in EEG, EOG, and EMG
- Amplitude may be increased due to impedance imbalance between inputs 1 and 2
- This will present often in obese patients because adipose tissue is much better conductor than muscle
 - Take care to not place M1, M2 on soft, fatty tissue
 - But still ECG artifact may not be avoidable
- Can link or “double reference” the 2 reference electrodes to minimize ECG artifact

ECG Artifact



Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.


Slow-Frequency Artifacts

- Usually result of perspiration or body movements associated with breathing
- Typically seen in EEG or EOG
- Very high amplitude, 1-2 Hz or slower in frequency
 - May look like N3 sleep
- Sweat artifact
 - Perspiration can loosen electrodes and dilute conductive gel/paste between electrode and skin
 - May need to reaffix electrodes and cool patient down
 - Cooling patient down may take several minutes
 - Salt bridge occurs when sweat communicates between 2 electrodes

Respiratory Artifact


- Typically recorded from electrodes around scalp and face
- Synchronous with patient's breathing rate and rhythm
- May be caused by:
 - Bed moving along with patient breathing
 - Unstable input board moving in sync with breathing
 - Electrode wires next to patient moving with each breath
 - Often patient is lying on the affected electrode
- May resolve with patient repositioning or replacing electrodes

Slow-Frequency Artifact



Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Slow-Frequency Artifact




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.


Pulse Artifact

- May occur in one or more EEG or EOG channels
- Occurs when electrode is placed over a pulsating blood vessel
- Will line up with QRS but produces sinusoidal slow wave rather than spike
- Need to relocate electrode

Pulse Artifact




From Butkov N: Atlas of clinical polysomnography, ed 2. Medford, Ore, 2010. Synapse Media. Reprinted with permission.



Cardioballistic Artifact

- Generated by cardiogenic oscillations, resulting in movement from intrathoracic pressure changes caused by contraction and relaxation of heart within thorax
- Not electrical artifact but lines up with QRS
- Typically on pressure and respiratory channels during central SDB event

Cardioballistic Artifact



Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Cardioballistic Artifact




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.


Channel-Blocking Artifact

- Occurs when amplitude of signal is higher than can be recorded within physical limitations of channel
- Often the result of faulty wire, popping electrode, or too high gain/sensitivity
- Creates a “flat” or “chopped off” appearance to the wave or signal may intrude into other channels
- If caused by profuse sweating, cool patient
 - If not caused by sweating, replace electrode or adjust gain/sensitivity
 - Last resort = Adjust LFF to 1 Hz

Channel-Blocking Artifact




Channel-Blocking Artifact



Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Channel-Blocking Artifact




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Electrode Popping Artifact


- Result of abrupt shift in electrode impedance
- Causes:
 - Mechanical instability of electrode
 - Intermittent short in electrode wire
 - Intermittently disconnected circuit
- Fixes:
 - Reaffix electrode
 - Secure circuit connection
 - Replace electrode
- LFF alters the duration of this type of artifact

Electrode Popping Artifact




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Electrode Popping Artifact



Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Electrode Popping Artifact




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Airflow Signal Artifacts

- Occurs when signal from thermistor or thermocouple degrades because of patient movement, sensor dislodged from path of airflow, oral breathing, or sensor contacting skin
- Must correct for malfunction as these channels imperative for respiratory event scoring and diagnosis of SDB

Airflow Signal Artifacts




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Respiratory Effort Sensor Artifact


- Causes:
 - Faulty sensor
 - Band displaced, over-tightened, or under-tightened
 - Connection unplugged or damaged
 - If external power, may be dead batteries
- Need to be efficient in isolating and correcting problem as these channels essential for scoring respiratory events and diagnosis SDB
 - Proper knowledge of equipment is essential

Respiratory Effort Sensor Artifact



Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Respiratory Effort Sensor Artifact




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Oximetry Artifacts


- Caused by faulty or absent pulse oximeter or bad connection
 - Also affects display of HR
- Oximetry essential for identification of hypopneas
- Other causes:
 - Poor perfusion
 - Severe hypoxemia
 - Nail polish or nail coverings
 - Movement

Oximetry Artifacts




Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

System Reference Artifact



From Butkov N. Atlas of clinical polysomnography, ed 2, Medford, Ore, 2010, Synapse Media. Reprinted with permission.

Disconnected Amplifier System/ Electrode Board



Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Disconnected Head Box



Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.

Implanted Electrical Devices

- Sometimes cranial implants and stimulators
 - May cause electrical interference
- Cardiac pacemakers
 - Ask patient if they have one
 - Small, battery-operated device typically implanted in chest just below right or left clavicle
- Atrial pacemaker
 - When SA node fails to consistently generate an impulse and circuitry to carry the signal is intact, patient only needs atrial pacemaker
 - Electrode rests in RA and stimulates atrium like SA node would

Atrial Pacemaker



From Aehlert B: ECGs made easy, ed 4, St Louis, 2011, Mosby JEMS-Elsevier.

Ventricular Pacemaker



From Aehlert B: ECGs made easy, ed 4, St Louis, 2011, Mosby JEMS-Elsevier.

AV Sequential Pacemaker



From Aehlert B: ECGs made easy, ed 4, St Louis, 2011, Mosby JEMS-Elsevier.

Fixed-Rate and Demand Pacemakers

- Fixed-rate generate impulse at specified rate and do not wait to sense if heart spontaneously generates own rhythm
- Demand sense whether heart is generating impulse and only fire when rate drops below programmed settings

Artifact Response

- Tech must determine when it is essential to disrupt patient's sleep to correct artifacts
 - Consider sleep stage, how long patient has been asleep, initial sleep latency, TST, time remaining, and whether there are back-up electrodes
 - Any time patient is awake, fix the artifact

Artifact Response

