

Surface Electrodes



- EEG typically amplified 20,000 times for visual display
- Electrodes = 1st and most important in patient-recorder interface
- Never mix electrodes made from different materials

Commonly Used Metals

- Most commonly used:
 - Gold-plated silver
 - Silver
 - Silver-silver chloride
- Used because of conductivity
- Application technique more important than the metal type
- Gold-plated electrodes are most commonly used in PSG
 - Should not be used if plating is chipped, cracked, or scratched because it causes artifact

Commonly Used Metals



Electrode Cups

- Stamped electrode cups are thinner and slightly bigger cup
 - Less durable
- Casted-style cups are thicker and smaller capacity
 - More durable
- Size
 - Adults = 10 mm in diameter
 - Pediatrics = 6 mm in diameter



Electrode Wires

- Also known as leads or lead wires
- Attach to cup electrodes
- Lead wires insulated with thin polytetrafluoroethylene are durable but tangle easy
- Thicker insulation may tangle less but not as sturdy



Specific Use for Electrodes



- Cup electrodes typically used for EEG and EOG
- Many labs use disposable stick-on patch electrodes for EMG and ECG
 - But can use cup electrodes for these parameters
- Facility decision to use disposable or reusable leads
 - Single use vs. multiple use

Electrode Connectors

- Must use sensors with a connector matching the sleep system used
- Connectors
 - Consult user manual for PSG system for proper connector
 - Most commonly used = 1.5-mm touchproof connectors
 - Also known as recessed-female connectors

Connectors



Used with permission from S&P
Whitaker Corporation, Inc., Sunnyvale,
Palo Alto, California.

1.5-mm
Touchproof



Used with permission from S&P
Whitaker Corporation, Inc., Sunnyvale,
Palo Alto, California.

2-pin 1-mm
Keyhole



Used with permission from S&P
Whitaker Corporation, Inc., Sunnyvale,
Palo Alto, California.

1/8-inch
(3.5-mm)
male phone

Connectors



RCA
connector



RJ11
connector



DIN
connector

Jumpers



Used with permission from MUSE Medical Supplies Inc, Menlo Park, California

- Typically used for “linked” mastoid references
- Allows signals to be referenced together to eliminate common artifact on both references
 - Helpful with ECG artifact
- Linking the mastoids also known as double referencing
 - Signals cancel each other out

Disposable Vs. Reusable Sensors

- Disposable
 - Single patient use
 - Do not require cleaning
 - More expensive
- Reusable
 - Multiple uses
 - Require cleaning and disinfection between applications
 - Less costly
 - Typically, higher quality signals



Cleaning Sensors

- Cleaning = Removal of foreign material from electrode cup
 - Uses water, detergent/soap, and mild scrubbing
 - First step before disinfection/sterilization
- Disinfection = Elimination of most pathogenic microorganisms (except bacterial spores)
 - Uses liquid chemicals (glutaraldehyde, sodium hypochlorite) or pasteurization (hot water)
- Sterilization = Complete elimination of all forms of microbial life, including spores
 - Uses physical or chemical processes such as autoclaving, ethylene oxide, or liquid chemicals

Cleaning Sensors

- Semicritical vs. Noncritical items
 - Noncritical = Electrodes and sensors that come in contact with intact skin and not mucous membranes
 - Require only low-level disinfection
 - Semicritical = Electrodes and sensors that come in contact with non-intact skin or mucous membranes
 - Require high-level disinfection

Interfacing Electrodes

- Derivations
 - Electrode inputs on head boxes form at least 2 types of derivations:
 - Referential = Measure difference in potential between 1 active electrode and 1 inactive electrode
 - Bipolar = Measure difference in potentials between 2 active electrodes
 - May be additional inputs for body temperature sensors, SpO2 sensors, nasal pressure devices, and other transducers

Operation, Use and Limitations of Surface Electrodes

- Site prep and proper electrode application essential for high-quality, low-noise recordings
- Improper cleaning of electrode sites leads to high impedances
 - Impedance = Combo of resistance and capacitance
 - Resistance = Electrical barrier
 - Capacitance = Ability of material to store electrical charge


Operation, Use and Limitations of Surface Electrodes

- Clean electrode site with alcohol to remove oil, then use abrading compound to remove dead skin cells
- Apply conductive paste/gel to electrode
- EEG = Recording of summed ionic flux from brain's neural activity
- EOG = Electrical changes occurring during eye movements as result of corneoretinal potential difference
- EMG = Recording of electrical impulses of muscle groups at rest and during contraction
- ECG = Electrical conduction through heart

Monitoring Airflow

- 2 most common methods to monitor airflow:
 - Thermal sensors
 - Detects airflow based on temperature changes
 - Used to detect apneas
 - Nasal pressure monitoring
 - Detects fluctuations in pressure caused by inspiration and expiration
 - Used to detect hypopneas and RERAs

Thermistor



Used with permission from MVAP Medical Supplies Inc. Newbury Park, California.

Monitoring Respiratory Effort

- Most common methods:
 - Respiratory inductive plethysmography (RIP)
 - Piezo technology
 - Diaphragmatic and intercostal EMG
- Other methods:
 - Mercury-filled strain gauges
 - Impedance pneumography


Piezo Technology


- Piezo-electric crystals are artificial or naturally occurring crystals that produce a charge output when compressed, flexed, or subjected to shear forces
 - When stress removed, electrical output stops
- Piezo belts measure tension where crystal is located
 - May have issues of accuracy when patient moves or false paradoxing



Diaphragmatic and Intercostal EMG

- Apply electrodes over intercostal spaces and upper abdomen
- Used to differentiate between central and noncentral respiratory events
- Reliability limited by proper electrode placement







Respiratory Inductive Plethysmography (RIP)

- Senses changes in cross-sectional area of rib cage and abdominal areas during breathing cycle
- Sensor consists of belt with wire woven or sewn in sine wave or zigzag pattern
- Creates magnetic field with module with circuit board, oscillator, and battery that passes weak current through wire
- More comfortable for patient
- Some include sum channel, useful for detecting paradoxical breathing


Respiratory Inductive Plethysmography (RIP)

- Signal can be calibrated or uncalibrated
 - Calibrated represent actual volume of airflow




Sum channel flattens or diminishes when chest and abdomen paradox.

RIP Vs. Piezo



Courtesy Mike Longman.

A close-up, black and white photograph of a microphone grille. The grille features a complex, woven mesh pattern with many small, rounded openings. The lighting creates highlights and shadows that emphasize the texture and depth of the mesh. The microphone is positioned on the left side of the frame, angled towards the right.

Snore Microphones and Sensors

- Microphone converts sound into small analog voltage
 - Most common types:
 - Electret
 - Dielectric material permanently electrically charged or polarized
 - Forms the diaphragm and its distance from plate causes voltage to be induced
 - Requires power source
 - Piezoelectric
 - Responds to vibrations on skin near upper airway
 - Does not require power source
 - Dynamic
 - Operate with moveable diaphragm displaced by sound wave
 - Movement creates voltage change proportional to power of sound wave
 - Does not require power source



Recording Body Position

- Important for determining if SDB is positional
- Can be monitored via video or through body position sensor
- Sensors provide data based on gravity
 - Must be placed in center of chest for accuracy

Audiovisual Monitoring and Recording

- Essential component of sleep study
- Video cameras are either fixed focus or pan-tilt-zoom
- Need infrared to see patient in dark
- 2 audio monitoring systems:
 - One-way communication from microphone to computer
 - On at all times to hear patient
 - Audio intercom system between control room and patient room

Differentiating Signal Types

- Bioelectric signals = Summed ionic flux generated by groups of cells or polarity of one anatomic location relative to another on an organ
 - Uses surface electrodes
- Transduced signals = Come from transducer that converts 1 type of energy to another
 - Converts mechanical energy to electrical energy
- Ancillary equipment = Any device that can process data on its own and interfaced with sleep system
 - EtCO₂ and SpO₂


Selecting an AC or DC Amplifier

- Rapidly fluctuating signals use AC amplifier
 - EEG, EOG, EMG, and ECG
 - Have HFF and LFF controls
- Slow-moving or constant variables use DC amplifier
 - PAP pressure readings and EtCO₂
 - Has HFF but not LFF

Sample Rate Requirements

- Analog-to-digital (A-to-D) conversion converts analog waveforms into series of numerical values for recording and display
- Sampling rate = Sampling frequency
 - Expressed in Hz
 - Higher sampling rate increases data resolution and file size

Signal Aliasing




Courtesy Mike Longman.

Original analog waveform

Sampled at 240



Sampled at 50




EEG Signals

- Require much higher sample frequency
- Nyquist-Shannon sampling theorem = Signal can be digitized and later restored to analog value if signal sampled at frequency greater than 2x the highest frequency contained in the signal
- AASM has specifications on minimal sample rates for PSG signals

EEG Signals



EEG Signals




Ancillary Recorder Signals, Equipment, and Methodology

- Esophageal manometry = Gold standard for assessing respiratory effort
 - Quantifies intrathoracic pressure changes
 - Not widely used in the clinical setting as it is invasive and may not be well tolerated
 - Very useful for detecting RERAs and diagnosing UARS

Recording and Understanding the Cyclic Alternating Pattern

- Also known as CAP
- Long-lasting periodic activity of 2 alternate EEG patterns consisting of transient arousals
- Divided into 2 phases:
 - Phase A
 - Greater arousal level
 - Phase B
 - Lesser arousal level
- Non-CAP = Prolonged stationary condition of sleep and arousal stability

CAP



A close-up image of a fingerprint, showing the intricate ridges and valleys. A red crosshair is overlaid on the image, consisting of a vertical line and a horizontal line intersecting in the center.

Overview of Pulse Oximetry

- Non-invasive and indirect method to monitor saturation of hemoglobin
- Expressed as SpO₂
- Sensor can be placed on fingertip, earlobe, toe, or forehead
 - Most sleep labs use fingertip sensors
- Shines light at 2 wavelengths
 - Red
 - Infrared

Overview of Pulse Oximetry

- Calculates light absorption at the 2 wavelengths to calculate portion of Hb saturated with O₂
- May cause inaccurate readings:
 - Poor perfusion
 - Anemia
 - CO poisoning
 - Movement
 - Nail polish
 - Dark-pigmented skin

Interfacing Ancillary Devices

- External device has output connection on back of device which gets connected to DC box connected to acquisition system
- Analog signal gets converted to digital signal
- Analog output = 0-1 V DC signal
- Each device has high and low output values
 - Low output = 0 V
 - High output = 1 V