



# Performing a CPAP Titration

---



# How CPAP Works

- Works by filtering normal room air (21% O<sub>2</sub>) and applying to a blower
- Measured in cm H<sub>2</sub>O



# Mask Fitting

- Many types and patient preference varies
- Review of types:
  - Nasal
  - Full Face
  - Pillows
  - Hybrid
  - Dental PAP
- Mask should be tight to skin but don't overtighten!

# CPAP Equipment

---

- Diagnostic units very different from home units
- Diagnostic units have different modes:
  - CPAP
  - Auto CPAP
  - Bi-level PAP
  - ASV
- Home units use interface device and store data on cards or on the Cloud
- Pressure relief
- Humidification

# Purpose of CPAP Titrations

- Titration = Gradual process of adjusting strength or dose of a medication or treatment until acceptable or optimal treatment level is achieved
- Patient hookup same as PSG except no thermistor or cannula
- Factors examined during titrations:
  - RDI
  - Mean and minimum SpO<sub>2</sub>
  - EEG arousal index
  - TST and REM duration at each treatment level

# Purpose of CPAP Titrations

- Goals:
  - RDI < 5
  - SpO2 > 90%
  - EEG arousals and RERAs < 5/hour
  - Snoring eliminated
- Goals may be difficult to achieve in single night
- Before titration, educate patient and try to desensitize them

# Purpose of CPAP Titrations

- PAP pressure verified by water column manometer
  - Air pushes water to certain level
- Patients may complain they don't sleep well the first night on CPAP
  - Document patient's subjective thoughts on CPAP
- Don't make treatment recommendations to patients
- Home setup for CPAP
  - Set up through DME company
  - AWAKE groups

# Increasing CPAP

- Follow lab protocols
- General rule:
  - Increase 2 cm at lower pressures
  - Increase 1 cm at higher pressures
  - Wait 5-20 min between pressure increases
- Increase for:
  - Obstructive apneas
  - Hypopneas
  - Snoring
  - EEG arousals
  - RERAs
  - O2 Desats



# Decreasing CPAP

- May need to decrease occasionally
- Typically at the beginning as patient tries to tolerate pressure
- Decrease for central apneas if higher pressure starts triggering them
- Document reason for decrease and when

# Bi-Level PAP

- Increases pressure during inhalation and decreases during exhalation
- IPAP = Pressure during inhalation
- EPAP = Pressure during exhalation
- Main purposes for use:
  - CSA
  - Patient intolerance of CPAP
  - Patient requires high level of PAP

# Bi-Level PAP

---

- CSA:
  - Apneas caused by hypocapnia (low CO<sub>2</sub>) or hypoxic respiratory drive combined with hyperoxemia (increased O<sub>2</sub> in blood)
- Bi-level PAP machines can become CPAPs by setting IPAP and EPAP at same level
- Settings:
  - Timed: Doesn't detect patient's changes in respirations
    - Reserved for severe patients
  - Spontaneous
  - Timed backup: Patient determines respiration changes until they don't breathe on own

# Bi-Level PAP

---

- Must keep a minimum of a 4 cm gap between IPAP and EPAP per AASM guidelines
- Purposes of IPAP and EPAP:
  - EPAP
    - Holds airway open
    - Increase to correct obstructive apneas
  - IPAP
    - Provides increased volume of air
    - Increase to increase baseline SpO2 and resolve flow limitation (hypopneas, snoring, RERAs, hypoventilation)

# Supplemental O2

---

May be used to increase baseline SpO2

---

Must have physician's order to add O2

---

May be used in conjunction with PAP

---

Usually used if low baseline O2 during wake

---

Titrate first for obstructions, then add O2 if necessary per lab protocol