

# Noninvasive Monitoring of Gas Exchange During Testing

Robertson Chapter 12

# Overview

- Involves obtaining O<sub>2</sub> or CO<sub>2</sub> values without skin puncture or invading the body
- Most common:
  - Pulse oximetry → O<sub>2</sub>Hb saturation
  - ETCO<sub>2</sub> or tcCo<sub>2</sub> for CO<sub>2</sub>
- Each device has limitations and advantages




# Overview of Gas Exchange

- Oxygenation values = Amount of O<sub>2</sub> carried by circulating blood
- Internal respiration = Process of gas diffusion between arterial blood and tissue at cellular level
- Diffusion = Movement of molecules, driven by a pressure gradient, across a semipermeable membrane

# Overview of Gas Exchange

---

- External respiration = Gases diffuse within the lungs across the alveolar-capillary membrane, where pulmonary capillaries and alveoli interface
- Inhaled air = 20.95% O<sub>2</sub>
- Respiration occurs only at cellular level
- Ventilation = Gross movement of air into and out of lungs
  - Moves gases between external environment and alveoli



# Hypoxia

---

Hypoxia = Inadequate tissue oxygenation

---

Hypoxemia = Low amounts of O<sub>2</sub> in blood

---

Hypoxic hypoxia = Tissue oxygenation diminished because of low amounts of O<sub>2</sub> in inspired air

---

Majority of O<sub>2</sub> transported by blood is carried to tissues in form of oxyhemoglobin (O<sub>2</sub>Hb)

---

Affinity of CO for Hb is about 210 times the affinity of O<sub>2</sub> for Hb

---

Anemic hypoxia = Hb impairment responsible for blood's reduced O<sub>2</sub> carrying capacity

# Hypoxia

- Circulatory hypoxia = Blood flow stagnant and doesn't deliver adequate O<sub>2</sub> to body
- Histotoxic hypoxia = Tissues become toxic and cells unable to use O<sub>2</sub> being delivered
  - Oxygenation appears better than actual results
- Lung disease can cause hypoxia
  - COPD
  - Obesity hypoventilation
  - Restrictive lung disorders

# Monitoring Oxygenation and Ventilation

---

- Oxygenation most commonly measured by ABG or pulse oximetry
  - Less commonly by transcutaneous
  - PaO<sub>2</sub> = ABG measurement
  - tcPO<sub>2</sub> = Transcutaneous measurement
  - SpO<sub>2</sub> = Pulse oximetry measurement
- Pulse oximeter = Non-invasive monitor for estimating saturation of Hb with O<sub>2</sub>
  - ABG value is calculated value not estimated like pulse oximetry
    - ABG = Gold standard of monitoring O<sub>2</sub> and CO<sub>2</sub>



# Oxygenation Status Classification

PaO <sub>2</sub>	SpO <sub>2</sub>	Oxygenation Status
80-100 mm Hg	95%-100%	Normal
60-80 mm Hg	90%-95%	Mild hypoxemia
40-60 mm Hg	70%-90%	Moderate hypoxemia
< 40 mm Hg	< 70%	Severe hypoxemia



# Monitoring CO<sub>2</sub>

- PaCO<sub>2</sub> = ABG measurement
- During hypoventilation, CO<sub>2</sub> in blood rises
  - Above 45 mm Hg = Hypercapnia (hypercarbia)
- During hyperventilation, CO<sub>2</sub> in blood decreases below 35 mm Hg
- 2 noninvasive measurements:
  - ETCO<sub>2</sub> = Sampling exhaled gas at end exhalation
  - tcPCO<sub>2</sub> = Transcutaneous

# Arterial Values of Blood Gases



- PaO<sub>2</sub> = Gold standard for determining oxygenation status
- Cooximetry = Additional test run in conjunction with ABG analysis that provides direct measurement of O<sub>2</sub>Hb saturation, abnormal Hb, and total Hb

Blood gas parameter	Arterial Value
pH	7.35-7.45
PaCO <sub>2</sub>	35-45 mm Hg
HCO <sub>3</sub> <sup>-</sup>	22-26 mEq/L
PaO <sub>2</sub>	80-100 mm Hg

# Assessment of pH Values

---

- Normal pH = 7.35-7.45
  - Even slight change can be life threatening
    - If  $< 6.8$  or  $> 7.8$ , cellular metabolism starts to stop, and death ensues
    - Metabolic component of acid-base balance assess based on  $\text{HCO}_3^-$ -value
      - Normally 22-26 mEq/L
      - $\text{HCO}_3^-$  changes show changes in renal system
  - Respiratory component is  $\text{PaCO}_2$ 
    - Hypoventilation =  $> 45$  mm Hg
    - Hyperventilation =  $< 35$  mm Hg

# Assessment of pH Values

---


- PaCo<sub>2</sub> = “Acid”
  - When increased, pH more acidic and decreases
  - When decreased, pH more alkaline and increases
- HCO<sub>3</sub><sup>-</sup> = “Base” (alkaline component)
  - When increased, causes upward shift in pH
  - When decreased, causes downward shift in pH
- Balance between acid and base responsible for pH value
  - Must be maintained for normal physiologic function

# Pulse Oximetry

---

- Often referred to as “fifth vital sign”
- Simple, noninvasive technique to monitor tissue oxygenation
- Data obtained through probe placed on each side of vascular bed (finger, ear, or toe)
  - Data represents % of saturated Hb with O<sub>2</sub>
- SpO<sub>2</sub> = Saturation of peripheral O<sub>2</sub>
  - Measured by pulse oximetry
- SaO<sub>2</sub> = Saturation of arterial O<sub>2</sub>
  - Measured by ABG

# Pulse Oximetry



*Modified from Gardner RM: J Cardiovasc Nurs 1:79, 1987.*

# Origins of Pulse Oximetry

---

- Dates back to early 1930s but no great interest until WWII
- Early 1940s – Glen Millikan developed lightweight ear O<sub>2</sub> meter for aviation monitoring
  - Coined term oximeter
  - Modified throughout 1940s and 1950s
  - Eventually manufactured by Waters Company
- 1970s – HP marketed system for use in clinical setting
  - Was bulky, heavy (35 lbs), and expensive (\$10,000)



# Origins of Pulse Oximetry

---

- 1972 – Modern oximeter developed by Takuo Aoyagi
- 1980 – Nellcor introduced first reliable and affordable oximeter commercially available
- Reusable sensors must be thoroughly cleaned and disinfected between patients
- There is a slight delay between actual oxygenation changes and the change in displayed value on pulse oximeter

# Reusable Adult Clip-Type Pulse Oximeter



From Elkin MK, et al: Nursing interventions and clinical skills, ed 4, St Louis, 2008, Elsevier.

# Pulse Oximeter Sensors



*From Linton AD: Introduction to medical-surgical nursing, ed 5, St Louis, 2012, Elsevier.*

# Theory of Pulse Oximetry

## Operation

---

- Use principles of spectrophotometry and photoplethysmography to estimate % of saturated Hb in blood
  - Spectrophotometry – Every substance has own pattern of light absorption
    - Pattern varies predictably with amount of substance that is present
    - Represents venous and capillary blood flow
  - Photoplethysmography – Uses light to detect minute volume changes between baseline and pulsatile component of blood flow through tissue bed
    - Represents arterial blood flow through tissue

# Theory of Pulse Oximetry

## Operation

---

- Uses 2 wavelengths of light: Red and infrared
- Transmitted light consists of 2 components:
  - DC component = Large, fixed component
    - Represents light passing through tissue and venous blood without being absorbed
  - AC component = Smaller portion pulsatile in nature
    - Changes absorption as blood pulses through arterioles
- Pulse oximeter divides AC signal by DC signal at each wavelength

# Oximetry Contraindications and Precautions

---

- Contraindications:
  - Presence of ongoing need for measurement of pH, PaCO<sub>2</sub>, Hbtot, and abnormal hemoglobins
  - Allergy to sensor material
- Precautions:
  - Burns due to prolonged use
  - Too tight → Pressure injury
  - Inspect device, power cord, and input cable for damage prior to use
  - Should be regularly tested for excessive electrical leakage current
  - Power cord should never be altered
  - Do not use incompatible probes

# Oximetry Infection Control

---

- Wear gloves
- Thoroughly clean and disinfect after each patient use for reusable probes
  - Also clean extension cable and oximeter itself after each use
- Throw away disposable probes after each use






# Limitations to Testing and Reliability of Results

---

- Techs should always attempt to correct problems and document what was tried
- Oxyhemoglobin Dissociation Curve
  - Oximeters estimate % of Hb that is saturated but cannot determine what Hb is saturated with
    - Can show a high O<sub>2</sub> reading even when patient has abnormally elevated levels of CO in Hb
    - May also present with clinical hypoxia even though showing normal % on oximeter
  - Peripheral perfusion, hypoxemia, hypotension, cold extremities, peripheral vasoconstriction, or reduction in diameter of blood vessels can make getting valid oximetry data difficult

# Plethysmographic Waveforms for Oximetry



*From Sandberg WS, et al: The MGH textbook of anesthetic equipment, ed 1, Philadelphia, 2011, Elsevier-Saunders.*


- A. Normal plethysmographic waveform
- B. Abnormal, unreliable plethysmographic waveform

# Additional Considerations for Pulse Oximetry Monitoring

---

- Nail coverings and nail polish should be removed prior to use
- Consider clinical issues that could potentiate increased motion artifact
- Make sure to properly place the sensor
  - Transducer and photodiode should lie directly across a vascular bed from one another
- Memory capacity, sampling rate, averaging time (response time), and data output capabilities can affect readings

# Factors Resulting in Erroneous Oximetry Values



Factor	Effect on SpO2 %
Nail polish	Falsely high
Dark skin pigmentation	Falsely high
Presence of COHb	Falsely high
High metHb	Falsely low if > 85% Falsely high if < 85%
Presence of fetal Hb	None
Anemia	None or low
Ambient light	Varies
Movement	Unpredictable
Poor perfusion	Poor signal; varies
Hypothermia	Poor signal; varies
Hypotension	Poor signal; varies

# Optimizing Oximetry During PSG

---

- AASM requires a max signal averaging time of 3 seconds for pulse oximetry during PSG
- Most oximeters interfaced with PSG acquisition systems provide 0-1 volt DC signal
- Tech must be proficient at troubleshooting equipment problems

# Recording Oximetry Data

---

- Operation manuals should be read for oximeters interfaced with PSG system
- Oximetry calibration
  - Simple procedure
  - Most oximeters generate calibration signals of 0 V and 1 V, representing 0% and 100% SpO<sub>2</sub>
- Selecting appropriate oximeter
  - Use the fastest response time setting available
  - CMS policy on supplemental O<sub>2</sub> is  $\leq 88\%$  for  $\geq 5$  minutes during sleep
  - Oximetry alone not sufficient for diagnosis
  - CMS and AASM require an associated O<sub>2</sub>Hb desat to score a hypopnea

# Monitoring Ventilation During Sleep

---

- Blood CO<sub>2</sub> levels provide most accurate representation of adequacy of ventilation
- To score hypoventilation in adults:
  - Increase in PaCO<sub>2</sub>  $\geq$  10 mm Hg during sleep to value exceeding 50 mm Hg for  $\geq$  10 minutes, or an increase to value  $>$  55 mm Hg for  $\geq$  10 minutes
  - Measured by ABG, ETCO<sub>2</sub>, or tcCO<sub>2</sub>
    - ETCO<sub>2</sub> and tcCO<sub>2</sub> primarily used in pediatrics



# Impact of Hypoventilation During Sleep

---

- Everyone has elevated CO<sub>2</sub> levels during sleep compared with wake values
- Normal adults experience increase in PaCO<sub>2</sub> of 2-8 mm Hg during NREM and 5-10 mm Hg during REM
- Underlying medical disease or SDB may cause patients to hypoventilate enough to affect pH
  - This compromises metabolic function
- Most labs don't monitor CO<sub>2</sub>
  - This can affect titration outcomes


# ETCO2 Monitoring

---


- Also known as capnometry or capnography
- Can provide info on patient's breathing pattern, underlying respiratory disease, and SDB events
- Often see increase in ETCO2 before desat in SpO2
- Often only used in pediatrics

# Monitoring CO<sub>2</sub> in Adults

- **Top:** Normal Pattern
- **Middle:** Hyperventilation
- **Bottom:** Ventilation-perfusion mismatching



From Ruppel GL: Manual of pulmonary function testing, ed 9, St Louis, 2009, Elsevier.




# Theory of Operation

---

- ETCO<sub>2</sub> monitor approximates partial pressure of CO<sub>2</sub> in arterial blood through end-tidal sampling of exhaled gas
- Infrared spectroscopy – Based on principle that molecules containing more than one element absorb infrared light in a characteristic manner
  - Can determine how much of a gas is present by measuring infrared radiation not absorbed



Microprocessor-Controlled Infrared CO2 Monitor



From Kacmarek RM, et al: Egan's fundamentals of respiratory care, ed 10, St Louis, 2013, Elsevier.

## Double-Beam Infrared Capnometer Schematic


# Theory of Operation

---

- 2 main methods of end-tidal monitoring:
  - Mainstream: Used when closed-system breathing circuit in use along with endotracheal or tracheostomy tube
  - Sidestream: For spontaneously breathing individuals
    - Typically, gas is sampled from nares through special nasal cannula
- Displaying ETCO<sub>2</sub> data
  - Single-breath and time-based capnography
    - Typically, time-based is used in PSG
  - Tech needs to make sure secretions, moisture, or organic material don't build up in tubing



# Capnogram



From Cairo JM: Pilbeam's mechanical ventilation: physiological and clinical applications, ed 5, St Louis, 2012, Elsevier.

# Advantages and Limitations

- Signal is not always reliable for sidestream sampling
- Mouth breathing can also be problematic
- Increased ventilation-to-perfusion mismatch will show CO<sub>2</sub> above the normal 4-6 mm Hg increase
  - Seen with severe COPD patients for example
- Some modern ETCO<sub>2</sub> monitors include integrated pulse oximetry
  - Good for showing trends

# Capnometer with Integrated Pulse Oximetry




From Cairo JM: Pilbeam's mechanical ventilation: physiological and clinical applications, ed 5, St Louis, 2012, Elsevier.

# Comparison of Obtained Values

---


- Good correlation between arterial, end-tidal, and transcutaneous CO<sub>2</sub> values
  - But ETCO<sub>2</sub> tends to underestimate CO<sub>2</sub> values obtained from arterial blood
    - Normal gradient is 4-6 mm Hg
    - Not a good tool when precise values are required
    - Diminished when there is increased flow from PAP or supplemental O<sub>2</sub>

# Hypnogram of Trend Capnography Data




From Beck SE, Marcus CL: Pediatric polysomnography, Sleep Med Clin 4:393-406, 2009.

# 7 Month Old with Meningomyelocele with Arnold-Chiari Malformation and Central Hypoventilation



From Beck SE, Marcus CL: Pediatric polysomnography. Sleep Med Clin 4:393-406, 2009.


# 4 Month Old Infant with OSA



From Beck SE, Marcus CL: Pediatric polysomnography. Sleep Med Clin 4:393-406, 2009.




# 6 Year Old with Obstructive Hypoventilation and Elevated ETCO2 with Prolonged Desat



From Beck SE, Marcus CL: Pediatric polysomnography. Sleep Med Clin 4:393-406, 2009.



# 60-Second Epoch of Periodic Breathing



From Beck SE, Marcus CL: Pediatric polysomnography, Sleep Med Clin 4:393-406, 2009.

# Contraindications, Precautions, and Infection Control

---

- Tech must carefully evaluate data obtained by sampling of end-tidal gas
- Gas sampling tubing and moisture filter are disposable, single-patient-use items
- No absolute contraindications for ETCO<sub>2</sub>


# Transcutaneous Monitoring of O<sub>2</sub> and CO<sub>2</sub>

---

- Approximates PaO<sub>2</sub> and PaCO<sub>2</sub> values by dilating capillary bed with heat and measuring capillary blood
- First mentioned by Severinghouse in 1960
  - First tcCO<sub>2</sub> sensor not commercially available until 20 years later
  - First tcO<sub>2</sub> sensor was used in 1967
  - Initially used in intensive care nursery
    - Sensors required frequent membrane changes and had to be moved to new site as often as every 2 hours
    - Monitor calibration required multiple times daily using gases of known low and high value

# Transcutaneous Monitoring of O<sub>2</sub> and CO<sub>2</sub>

- Typically only used in pediatrics
- Today's sensors can remain in one location as long as 8-12 hours
  - Sensors are small and lightweight



From Cairo JM, Pilbeam SP: Mosby's respiratory care equipment, ed 9, St Louis, 2011, Elsevier.

# Advantages and Limitations

---

- Measured on surface of skin using modified blood gas electrodes
- Optimal location for sensor placement is area of high capillary density and minimal skin thickness
  - Earlobe, forearm, chest and abdomen
- Advantages:
  - No need to draw arterial blood
  - Ability to monitor and obtain reliable results during PAP or supplemental O<sub>2</sub>
- Disadvantage: Need for recording site to be warmed by sensor in order to enhance capillary vasodilation for increased accuracy
  - Can take 12-30 minutes for site to warm up

# Comparison of Obtained Values


- tcCO<sub>2</sub> monitoring more reliable than tcO<sub>2</sub> monitoring
  - tcO<sub>2</sub> provides little value in a sleep study
- tcCO<sub>2</sub> tends to slightly overestimate PaCO<sub>2</sub>
- tcCO<sub>2</sub> provides most reliable data for determining effects of therapy on ventilation throughout sleep study
  - Provides valid results when ETCO<sub>2</sub> may not be reliable

# Contraindications, Precautions, and Infection Control

---

- No known absolute contraindications to tcCO<sub>2</sub> monitoring
- Must be careful not to leave probe in one spot too long so as not to damage skin
- Integrity of skin must be assessed prior to placement of sensor
- Sensor, cable, and recorder should be cleaned and disinfected between each patient

# Acid-Base with Normal pH




\* The other parameter will be moving in the opposite direction of the pH (acid or alkaline) to compensate

Copyright © 2014 by Mosby, Inc., an affiliate of Elsevier Inc.



# Acid-Base with Abnormal pH



# Acid-Base Balance Quick Reference

	Acute	Partially Compensated	Fully Compensated
<b>Respiratory Acidosis</b>			
pH	↓	↓	Returns to normal
PaCO <sub>2</sub>	↑	↑	↑
HCO <sub>3</sub> <sup>-</sup>	Normal	↑	↑
<b>Respiratory Alkalosis</b>			
pH	↑	↑	Returns to normal
PaCO <sub>2</sub>	↓	↓	↓
HCO <sub>3</sub> <sup>-</sup>	Normal	↓	↓

# Acid-Base Balance Quick Reference

	Acute	Partially Compensated	Fully Compensated
<b>Metabolic Acidosis</b>			
pH	↓	↓	Returns to normal
PaCO <sub>2</sub>	Normal	↓	↓
HCO <sub>3</sub> <sup>-</sup>	↓	↓	↓
<b>Metabolic Alkalosis</b>			
pH	↑	↑	Returns to normal
PaCO <sub>2</sub>	Normal	↑	↑
HCO <sub>3</sub> <sup>-</sup>	↑	↑	↑

# Acid-Base Balance Quick Reference

	Combined Acidosis	Combined Alkalosis
pH	↓	↓
PaCO <sub>2</sub>	↑	↓
HCO <sub>3</sub> <sup>-</sup>	↓	↑